

# Comfort-fit Face Mask



## The Next Generation in Comfort

### ASTM F2100-23 Level 1 / 2 / 3

#### Product Features

- Comfort-fit with 3-ply medical-grade protection
- Neatly prefolded, anti-fatigue ribbon ear-loops
- Dual core moldable nose-piece for secure fit
- ASTM F2100 Level 1 / 2 / 3 certifications



**Front**



**Back**



Anti-Fatigue Ear-loops

**MADE IN CANADA**



**Canadian-made at Competitive Prices**

- Manufactured in Ontario
- Stock available for quick delivery
- Competitive pricing



**NO LATEX / NON-STERILE / SINGLE USE**

#### Specifications

Material	Mask Body	Outer Layer	Spunbond Polypropylene (SBPP)
		Filtration Layer	Meltblown Polypropylene (MBPP)
		Inner Layer	Spunbond Polypropylene (SBPP)
Dimensions	Nose-piece	Aluminium	
	Ear-loops	Ethylene-based polymer & Polypropylene	
	Body Size	Adult - 175mm x 95mm	Child - 145mm x 95mm

#### ASTM F2100-23 Requirements

	LEVEL 1	LEVEL 2	LEVEL 3
Bacterial filtration efficiency (BFE)	≥95%	≥98%	≥98%
Sub-micron particulate filtration efficiency at 0.1 µm (PFE)	≥80%	≥85%	≥85%
Resistance to penetration by synthetic blood [mm Hg]	@80	@120	@160
Differential Pressure (ΔP) [mm H <sub>2</sub> O/cm <sup>2</sup> ]	< 5.0	< 6.0	< 6.0
Flame Spread	Class 1	Class 1	Class 1



[www.vmedcare.co](http://www.vmedcare.co)

VMedCare™ is the PPE division of Viva Healthcare Packaging



[www.viva-healthcare.com](http://www.viva-healthcare.com)