

PPE TURNKEY SOLUTIONS

Personal Protective Equipment



- Surgical Face Masks
- Face Shields
- Hand Sanitizers
- Alcohol Wipes
- Disposable Gloves
- Custom Kit Options



MADE IN CANADA

PPE items manufactured by Viva in Canada



Viva Healthcare Packaging

www.vmedcare.co





Surgical Face Masks (3-layer, pleated, ear-loop)

- Dimensions: Adult 175 x 95mm, Child 145 x 95mm
- ASTM F2100 Level 1 / Level 2 / Level 3

DESCRIPTION	PART NO.	PACKING
ASTM F2100 Level 1 Face Mask For Adult	VM-AL01	25pcs / 50pcs per box
ASTM F2100 Level 2 Face Mask For Adult	VM-AL02	25pcs / 50pcs per box
ASTM F2100 Level 3 Face Mask For Adult	VM-AL03	25pcs / 50pcs per box
ASTM F2100 Level 1 Face Mask For Child	VM-AL01K	50pcs per box
ASTM F2100 Level 2 Face Mask For Child	VM-AL02K	50pcs per box
ASTM F2100 Level 3 Face Mask For Child	VM-AL03K	50pcs per box



Face Shields

- Medical Face Shield (disposable) – maximum safety in medical environment
- Retractable Face Shield (reusable) – visor can be lifted up 90° for non-medical environment
- Custom pack options for single visor and multiple replacement lenses
- Anti-fog and anti-static lens treatment

DESCRIPTION	PART NO.	PACKING
Medical Face Shield	L600	100 pcs per box
Retractable Face Shield (10 Lens)	VM-RE	30 pcs per box



Hand Sanitizers

- Active ingredient = 65% Ethyl Alcohol
- Sachet size: 60 x 70mm

DESCRIPTION	PART NO.	PACKING
Antibacterial Hand Sanitizer 3mL	VM-HS3	50/100/1,000 sachets per box



Alcohol Wipes

- Saturated with 64% Isopropyl Alcohol
- Non-woven cloth
- Sachet size: 100 x 60mm
- Wipe size: 180 x 120mm

DESCRIPTION	PART NO.	PACKING
64% Alcohol Wipe	VM-AW	50/100/1,000 sachets per box



Disposable Gloves

- Available upon request
- Viva to outsource from qualified manufacturer

DESCRIPTION	PART NO.	PACKING
Nitrile Glove (L size)	VNGLOVE-L	100pcs per box
Nitrile Glove (XL size)	VNGLOVE-XL	100pcs per box
PVC Glove (L size)	VPGLOVE-L	100pcs per box
PVC Glove (XL size)	VPGLOVE-XL	100pcs per box



Custom Kit Options

- Traveller kits in custom-printed polybag with zipper lock
- Box set in carton or rigid plastic box for more PPE items
- Kits assembled by Viva in Canada