

Gen2 Face Mask


The Next Generation in Comfort

ASTM F2100 Level 1 / 2 / 3



Product Features

- Comfort fit with 3-ply protection
- Neatly prefolded, easy-tear and stretchable ear-loops
- Dual core moldable nose-piece for secure fit
- ASTM F2100 Level 1 / 2 / 3 certifications



Outside View



Inside View

MADE IN CANADA



Canadian-made at Competitive Prices

- Manufactured in Ontario
- Stock available for quick delivery
- Health Canada approved
- Competitive pricing



NO LATEX / NON-STERILE / SINGLE USE

Materials

Mask Body	Outer Layer	Spunbond Polypropylene (SBPP)
	Filtration Layer	Meltblown Polypropylene (MBPP)
	Inner Layer	Spunbond Polypropylene (SBPP)
Nose-piece	Aluminium	
Ear-loop	Ethylene-based polymer & Polypropylene	

ASTM F2100 Requirements

	LEVEL 1	LEVEL 2	LEVEL 3
Bacterial filtration efficiency (BFE)	≥95%	≥98%	≥98%
Sub-micron particulate filtration efficiency at 0.1 µm (PFE)	≥95%	≥98%	≥98%
Resistance to penetration by synthetic blood [mm Hg]	@80	@120	@160
Differential Pressure (ΔP) [mm H ₂ O/cm ²]	< 5.0	< 6.0	< 6.0
Flame Spread	Class 1	Class 1	Class 1



www.vmedcare.co

VMedCare™ is the PPE division of Viva Healthcare Packaging



www.viva-healthcare.com